

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 1 of 8

---

## Section 1 - CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

---

### PRODUCT NAME

DENTSPLY DIE COLOUR LAYER ISOLATION

### SYNONYMS

! 01/02

### PRODUCT USE

Dental restoration.

### SUPPLIER

Company: Dentsply (Australia) Pty Ltd (ABN: 15 004 290 322)

Address:  
11-21 Gilby Road  
Mount Waverley  
VIC, 3149  
AUS

Telephone: +61 03 9538 8240

Emergency Tel: 0413 830 239

Fax: 03 9538 8260

---

## Section 2 - HAZARDS IDENTIFICATION

---

### STATEMENT OF HAZARDOUS NATURE

NON-HAZARDOUS SUBSTANCE. NON-DANGEROUS GOODS.

According to the Criteria of NOHSC, and the ADG Code.

### POISONS SCHEDULE

None

### RISK

### SAFETY

---

## Section 3 - COMPOSITION / INFORMATION ON INGREDIENTS

---

NAME	CAS RN	%
surfactant		NotSpec
water	7732-18-5	NotSpec

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 2 of 8

---

## Section 4 - FIRST AID MEASURES

---

### SWALLOWED

- Immediately give a glass of water.
- First aid is not generally required. If in doubt, contact a Poisons Information Centre or a doctor.

### EYE

- If in eyes, hold eyelids apart and flush the eye continuously with running water.
- Continue flushing until advised to stop by the Poisons Information Centre or a doctor, or for at least 15 minutes.
- Ensure complete irrigation of the eye by keeping eyelids apart and away from eye and moving the eyelids by occasionally lifting the upper and lower lids.
- If pain persists or recurs seek medical attention.
- Removal of contact lenses after an eye injury should only be undertaken by skilled personnel.

### SKIN

If skin or hair contact occurs:

- Flush skin and hair with running water (and soap if available).
- Seek medical attention in event of irritation.

### INHALED

- If fumes or combustion products are inhaled remove from contaminated area.
- Other measures are usually unnecessary.

### NOTES TO PHYSICIAN

Treat symptomatically.

---

## Section 5 - FIRE FIGHTING MEASURES

---

### EXTINGUISHING MEDIA

- Water spray or fog.
- Foam.
- Dry chemical powder.
- BCF (where regulations permit).
- Carbon dioxide.

### FIRE FIGHTING

- Alert Fire Brigade and tell them location and nature of hazard.
- Wear full body protective clothing with breathing apparatus.
- Prevent, by any means available, spillage from entering drains or water course.
- Use water delivered as a fine spray to control fire and cool adjacent area.
- Avoid spraying water onto liquid pools.

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 3 of 8

---

## Section 5 - FIRE FIGHTING MEASURES ...

---

- DO NOT approach containers suspected to be hot.
- Cool fire exposed containers with water spray from a protected location.
- If safe to do so, remove containers from path of fire.

### FIRE/EXPLOSION HAZARD

- Combustible.
  - Slight fire hazard when exposed to heat or flame.
  - Heating may cause expansion or decomposition leading to violent rupture of containers.
  - On combustion, may emit toxic fumes of carbon monoxide (CO).
  - May emit acrid smoke.
  - Mists containing combustible materials may be explosive.
- Combustion products include  
carbon dioxide (CO<sub>2</sub>)  
other pyrolysis products typical of burning organic material

### FIRE INCOMPATIBILITY

Avoid contamination with oxidising agents i.e. nitrates, oxidising acids, chlorine bleaches, pool chlorine etc. as ignition may result

### HAZCHEM

None

---

## Section 6 - ACCIDENTAL RELEASE MEASURES

---

### EMERGENCY PROCEDURES

#### MINOR SPILLS

- Remove all ignition sources.
- Clean up all spills immediately.
- Avoid breathing vapours and contact with skin and eyes.
- Control personal contact by using protective equipment.
- Contain and absorb spill with sand, earth, inert material or vermiculite.
- Wipe up.
- Place in a suitable labelled container for waste disposal.

#### MAJOR SPILLS

Moderate hazard.

- Clear area of personnel and move upwind.
- Alert Fire Brigade and tell them location and nature of hazard.
- Wear breathing apparatus plus protective gloves.
- Prevent, by any means available, spillage from entering drains or water course.
- No smoking, naked lights or ignition sources.
- Increase ventilation.
- Stop leak if safe to do so.
- Contain spill with sand, earth or vermiculite.

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 4 of 8

---

## Section 6 - ACCIDENTAL RELEASE MEASURES ...

---

- Collect recoverable product into labelled containers for recycling.
- Absorb remaining product with sand, earth or vermiculite.
- Collect solid residues and seal in labelled drums for disposal.
- Wash area and prevent runoff into drains.
- If contamination of drains or waterways occurs, advise emergency services.

**Personal Protective Equipment advice is contained in Section 8 of the MSDS.**

---

## Section 7 - HANDLING AND STORAGE

---

### PROCEDURE FOR HANDLING

- Avoid all personal contact, including inhalation.
- Wear protective clothing when risk of exposure occurs.
- Use in a well-ventilated area.
- Prevent concentration in hollows and sumps.
- DO NOT enter confined spaces until atmosphere has been checked.
- Avoid smoking, naked lights or ignition sources.
- Avoid contact with incompatible materials.
- When handling, DO NOT eat, drink or smoke.
- Keep containers securely sealed when not in use.
- Avoid physical damage to containers.
- Always wash hands with soap and water after handling.
- Work clothes should be laundered separately.
- Use good occupational work practice.
- Observe manufacturer's storing and handling recommendations.
- Atmosphere should be regularly checked against established exposure standards to ensure safe working conditions.

### SUITABLE CONTAINER

- Metal can or drum
- Packaging as recommended by manufacturer.
- Check all containers are clearly labelled and free from leaks.

### STORAGE INCOMPATIBILITY

Avoid contamination of water, foodstuffs, feed or seed.  
Avoid reaction with oxidising agents

### STORAGE REQUIREMENTS

- Store in original containers.
- Keep containers securely sealed.
- Store in a cool, dry, well ventilated area.
- DO NOT allow to freeze.
- Store away from incompatible materials.
- Protect containers against physical damage and check regularly for leaks.
- Observe manufacturer's storing and handling recommendations.

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 5 of 8

---

## Section 8 - EXPOSURE CONTROLS / PERSONAL PROTECTION

---

### EXPOSURE CONTROLS

None assigned. Refer to individual constituents.

### INGREDIENT DATA

#### WATER:

No exposure limits set by NOHSC or ACGIH

### PERSONAL PROTECTION

#### EYE

- Safety glasses with side shields
- Chemical goggles.
- Contact lenses pose a special hazard; soft lenses may absorb irritants and all lenses concentrate them.

#### HANDS/FEET

Wear general protective gloves, eg. light weight rubber gloves.

#### OTHER

No special equipment needed when handling small quantities.

#### OTHERWISE:

- Overalls.
- Barrier cream.
- Eyewash unit.

The local concentration of material, quantity and conditions of use determine the type of personal protective equipment required. For further information consult site specific CHEMWATCH data (if available), or your Occupational Health and Safety Advisor.

### ENGINEERING CONTROLS

General exhaust is adequate under normal operating conditions. If risk of overexposure exists, wear SAA approved respirator. Correct fit is essential to obtain adequate protection. Provide adequate ventilation in warehouse or closed storage areas.

---

## Section 9 - PHYSICAL AND CHEMICAL PROPERTIES

---

### APPEARANCE

Colourless liquid; mixes with water.

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 6 of 8

---

## Section 9 - PHYSICAL AND CHEMICAL PROPERTIES ...

---

### PHYSICAL PROPERTIES

Liquid.  
Mixes with water.

Molecular Weight: Not Applicable  
Melting Range (°C): Not Available  
Solubility in water (g/L): Miscible  
pH (1% solution): Not Available  
Volatile Component (%vol): Not Available  
Relative Vapour Density (air=1): Not Available  
Lower Explosive Limit (%): Not Available  
Autoignition Temp (°C): Not Available  
State: Liquid

Boiling Range (°C): Not Available  
Specific Gravity (water=1): 0.96 @ 20 degC  
pH (as supplied): 4-5  
Vapour Pressure (kPa): Not Available  
Evaporation Rate: Not Available  
Flash Point (°C): Not Available  
Upper Explosive Limit (%): Not Available  
Decomposition Temp (°C): Not Available

---

## Section 10 - CHEMICAL STABILITY AND REACTIVITY INFORMATION

---

### CONDITIONS CONTRIBUTING TO INSTABILITY

- Presence of incompatible materials.
- Product is considered stable.
- Hazardous polymerisation will not occur.

---

## Section 11 - TOXICOLOGICAL INFORMATION

---

### POTENTIAL HEALTH EFFECTS

#### ACUTE HEALTH EFFECTS

#### SWALLOWED

(No Oral LD50, any animal species) The material has NOT been classified by EC Directives or other classification systems as "harmful by ingestion". This is because of the lack of corroborating animal or human evidence. The material may still be damaging to the health of the individual, following ingestion, especially where pre-existing organ (eg. liver, kidney) damage is evident. Present definitions of harmful or toxic substances are generally based on doses producing mortality rather than those producing morbidity (disease, ill-health). Gastrointestinal tract discomfort may produce nausea and vomiting. In an occupational setting however, ingestion of insignificant quantities is not thought to be cause for concern.

#### EYE

Although the liquid is not thought to be an irritant (as classified by EC Directives), direct contact with the eye may produce transient discomfort characterised by tearing or conjunctival redness (as with windburn).

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 7 of 8

---

## Section 11 - TOXICOLOGICAL INFORMATION ...

---

### SKIN

The material is not thought to produce adverse health effects or skin irritation following contact (as classified by EC Directives using animal models). Nevertheless, good hygiene practice requires that exposure be kept to a minimum and that suitable gloves be used in an occupational setting.

### INHALED

The material is not thought to produce adverse health effects or irritation of the respiratory tract (as classified by EC Directives using animal models). Nevertheless, good hygiene practice requires that exposure be kept to a minimum and that suitable control measures be used in an occupational setting.

### CHRONIC HEALTH EFFECTS

Long-term exposure to the product is not thought to produce chronic effects adverse to the health (as classified by EC Directives using animal models); nevertheless exposure by all routes should be minimised as a matter of course.

## Dentsply Die Colour Layer Isolation

Not available. Refer to individual constituents.  
unless otherwise specified data extracted from RTECS - Register of Toxic Effects of Chemical Substances

### WATER:

No significant acute toxicological data identified in literature search.

---

## Section 12 - ECOLOGICAL INFORMATION

---

---

## Section 13 - DISPOSAL CONSIDERATIONS

---

- Recycle wherever possible or consult manufacturer for recycling options.
- Consult State Land Waste Management Authority for disposal.
- Bury residue in an authorised landfill.
- Recycle containers if possible, or dispose of in an authorised landfill.

---

## Section 14 - TRANSPORTATION INFORMATION

---

Shipping Name:

None

Dangerous Goods Class: None

UN/NA Number: None

ADR Number: None

Packing Group: None

Labels Required:

Additional Shipping Information:

International Transport Regulations:

IMO: None

continued...

# DENTSPLY DIE COLOUR LAYER ISOLATION

ChemWatch Material Safety Data Sheet  
Issue Date: Thu 14-Oct-2004

CHEMWATCH 4613-97  
CD 2004/3 Page 8 of 8

---

## Section 14 - TRANSPORTATION INFORMATION ...

---

### HAZCHEM

None

---

## Section 15 - REGULATORY INFORMATION

---

### POISONS SCHEDULE

None

### REGULATIONS

The following substances are found on Australian Inventory of Chemical Substances (AICS):  
water (CAS: 7732-18-5)

---

## Section 16 - OTHER INFORMATION

---

This document is copyright. Apart from any fair dealing for the purposes of private study, research, review or criticism, as permitted under the Copyright Act, no part may be reproduced by any process without written permission from CHEMWATCH. TEL (+61 3) 9572 4700.

Issue Date: Thu 14-Oct-2004  
Print Date: Wed 20-Oct-2004